



**URGENT MEDICAL DEVICE  
RECALL**

**MAKING LIFE EASY FOREARM  
CRUTCH (TALL / ADULT)**

**Batch numbers:**

**MLE Forearm Crutch Tall: 251097-2105-1, 251097-2107-1,  
251097-2201-1**

**MLE Forearm Crutch Adult: 234354-2105-1, 234354-2107-1,  
234354-2201-1, 234354-2204-1**

Australian Pharmaceutical Pty Ltd, following consultation with the Therapeutic Goods Administration, is recalling the above stated batches of MAKING LIFE EASY Forearm Crutch (Tall / Adult (which is a walking aid) as the hand grip on the handle may detach while using the product. No other batches of Forearm crutch are affected by this recall.

If you have a MAKING LIFE EASY Forearm crutch from the stated batches, do not use it. Return it to the place of purchase for a refund or call customer service line to arrange the return of affected product and refund.

**CUSTOMER SERVICE 1300 363 303**

Using the defective product may result in a fall causing injury. Anyone who is concerned in any way about the use of this product should consult their service provider.

Australian Pharmaceutical Pty Ltd sincerely regrets any inconvenience to their customers.